ceramill® matik

ceramill® cleanstream

DE	Betriebsanleitung	3 - 20
EN	User Manual	21 - 37
FR	Mode d'emploi	38 - 56
IT	Istruzioni d'uso	57 - 75
ES	Modo de empleo	76 - 94
PT	Instruções de operação	95 - 113

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Tab. 2

- Translation of the original Instruction Manual -

Table of Contents

1	Explanation of Symbols	22	9	Malfunctions, Repairs and Warranty ..	35
			9.1	Malfunctions	35
			9.2	Repairs	35
			9.3	Warranty	35
2	Safety Instructions	22	10	Environmental Protection	35
2.1	General Safety Instructions	22	10.1	Packaging	35
2.2	Unsupervised operation	23	10.2	Waste equipment	35
2.3	Machining of titanium blanks	24	10.3	Coolant/Lubricant (C/L)	35
3	Suitable Personnel	24	11	Technical Data	36
4	Information on the System	24			
4.1	General	24			
4.2	Delivery Scope	25			
4.3	Intended Use	25			
4.4	CE Declaration of Conformity	26			
4.5	Components and interfaces	27			
4.6	Accessories	30			
5	Installation	31			
5.1	Country-specific supply voltage ...	31			
5.2	Setting up	31			
5.3	Maintenance unit and air purity during operation	32			
5.4	Requirements for the extraction system	32			
6	Application and Operation	32			
7	Handling the Coolant/Lubricant (C/L) .	33			
8	Cleaning and Maintenance	34			

1 Explanation of Symbols

Warning indications

Warning indications in the text are marked with a triangle and boxed.

In case of hazards through electricity, the exclamation mark in the warning triangle is substituted by a lightning bolt.

Signal words at the beginning of a warning indication specify the type and severity of the consequences, if the measures to avert the hazard are not adhered to.

- _ **NOTE** means that property damage can occur.
- _ **CAUTION** means that light to fairly serious personal injury can occur.
- _ **WARNING** means that serious personal injury can occur.
- _ **DANGER** means that life-threatening personal injury can occur.

Important information

Important information that do not lead to hazards for humans or property damage are marked with the icon aside and are also boxed.

Other symbols in the Manual

Symbol	Meaning
▷	Item of an operation description
—	Item of a list
•	Subitem of an operation description or a list
[3]	Numbers in square brackets refer to position numbers in graphics/figures

2 Safety Instructions

2.1 General Safety Instructions

When using the product, always observe the following safety instructions:

DANGER:

Danger of electric shock from leaking coolant/lubricant.

- ▷ When changing the coolant/lubricant, use the intended pumping device.
- ▷ When coolant/lubricant has leaked: Disconnect the machine from the power supply and contact the technical service.

CAUTION:

Damage of slipping from leaking coolant/lubricant!

- ▷ When changing the coolant/lubricant, use the intended pumping device.
- ▷ When coolant/lubricant has leaked: Set up warning signs and clean floor properly.

CAUTION:

Danger of tipping over when the drawer for the belt filter and tank is pulled out.

- ▷ Do not use the drawer as a climbing aid or seat.

CAUTION:

A defective product can lead to malfunctions!

On detection of damage or a functional defect of the product:

- ▷ Label the product as defective.
- ▷ Prevent further operation until the product has been repaired.

**NOTE:**

The system may only be used for the purpose that it is intended for. When disregarding, damage to property and personal injury cannot be excluded. Any and all warranty, guarantee and responsibility on behalf of the manufacturer shall expire.

**NOTE:**

The system may only be operated in enclosed rooms.

**NOTE:**

Possible damage due to insufficient load-bearing capacity of the floor.

The set-up location and the surrounding paths must be able to withstand the weight of the system.

- ▷ Observe weight of system.
- ▷ Inform yourself about the building conditions at the responsible authority (e.g. architect/structural engineer).

**NOTE:**

Possible machine damage through insufficient extraction during dry processing!

- ▷ Operate the machine only with the extraction system recommended by the manufacturer or another compatible extraction system.

**NOTE:**

Possible machine damage due to insufficient or incorrect compressed-air supply!

- ▷ Operate the machine with the required pressure, air volume and air quality (see Technical Data).

**NOTE:**

Switch the machine off when not in use or unsupervised for longer periods, e.g., overnight. This method of saving electrical energy also benefits the environment.

2.2 Unsupervised operation

**NOTE:**

The machine may be operated with no supervision, provided that the national and local laws and regulations permit this and that these are complied with. Exception: Machining of titanium blanks (see section 2.3). Furthermore, the requirements of the respective insurance company must be met.

The following conditions must be met to ensure unsupervised operation:

- ▷ Do not use the machine with a heavily contaminated milling chamber.
- ▷ Protect the machine against interference by unauthorized persons.
- ▷ Ensure that the room in which the machine is being operated has an automatic fire alarm system.

2.3 Machining of titanium blanks

NOTE:

Machining titanium involves an increased risk of fire despite water cooling. Metal chips can easily ignite and cause a metal fire.

For titanium indications:

- ▷ Do not operate machine unsupervised.
- ▷ Have suitable extinguishing means readily available.

3 Suitable Personnel

NOTE:

The machine may only be unpacked, installed or relocated by trained specialist personnel, authorized by the manufacturer.

NOTE:

The machine may only be put into operation by trained specialist personnel, authorized by the manufacturer.

NOTE:

The machine may only be operated by trained specialised personnel.

4 Information on the System

4.1 General

This User Manual describes various models of the system. Therefore, the representations can deviate from the actual appearance of the machine.

The system is provided factory-calibrated. Prior to delivery, a test object was factory-milled in order to check the zero point settings. Possibly remaining milled shavings/chips are not a reason for rejection.

- ▷ After unpacking by the authorised shipping partner: check for possible transport damage. In case of transport damage, immediately file complaint on site with the freight carrier.
- ▷ After installation by the authorised service technician: Check scope of delivery and correct function. In case of missing parts or incorrect function, immediately file complaint on site with the service technician.

4.2 Delivery Scope

Ceramill Matik:

- _ Milling machine (Ceramill Matik)
- _ Key for manual unlocking of the milling compartment
- _ Scanner Voyager, set incl. operating instructions and USB connection cable
- _ Service set, spindle
- _ Maintenance unit incl. connection material
- _ 2 x Torque wrench 2 Nm
- _ Interchangeable blade, size 3
- _ Interchangeable blade, size 2.5
- _ Allen key with T-handle, size 3
- _ Angle Allen key, size 2.5
- _ Angle Allen key, size 5
- _ Calipers 0 - 150 mm
- _ 2 test blanks for calibration milling
- _ Wax blank
- _ Cleaning brush (large)
- _ Cleaning brush (small)
- _ Control line, extraction
- _ USB cable 3 m
- _ Network cable 5 m
- _ Power cable
- _ Workpiece holder 98 (clamped in machine)
- _ Tool starter set
- _ Glass-tube fuse set
- _ Connection set for service technician
- _ Adapter for extraction
- _ Suction strainer

Ceramill Cleanstream:

- _ Coolant/lubricant pumping system with integrated chip-removal system (Ceramill Cleanstream)
- _ Filter fleece for belt filter incl. roller unit
- _ Filter cartridge 10" 5µm
- _ Analytics instructions
- _ KSS test record
- _ 3 x 1 l "Coolant/lubricant"
- _ pH test-strip set
- _ Nitrate test set

- _ Refractometer for determination of coolant content
- _ Towel
- _ Filter spanner, size 2
- _ Allen key, size 10
- _ Power cable
- _ USB cable 30 cm
- _ Cleaning gun with spiral hose
- _ Piping kit

4.3 Intended Use

The system is a PC-controlled milling and grinding machine for the fabrication of dentures by means of:

Ceramill Matik

- _ Dry processing of blanks made of
 - Pre-sintered zirconium oxide
 - Unsintered non-precious metals (NPMs)
 - Wax
 - Plastics
 - Hybrid materials
 - Modelling materials

Ceramill Matik in conjunction with Ceramill Cleanstream

- _ Dry processing of blanks made of
 - Cobalt-chromium alloys
- _ Wet processing of blanks made of
 - Sintered glass ceramics
 - Plastics
 - Hybrid materials
 - Feldspathic ceramics
 - Titanium
 - Cobalt-chromium alloys

Using blanks and tools not approved by the manufacturer can damage the machine and make the product unusable. For such cases, the manufacturer shall assume no liability whatsoever. Unauthorised modifications/alterations of the machine shall also void the warranty.

Operation of the system is permitted only with a suitable extraction system. Consult the manufacturer for more information (for technical requirements see chapter 5.4).

NOTE:

The system may only be used for the purpose that it is intended for. When disregarding, damage to property and personal injury cannot be excluded. Any and all warranty, guarantee and responsibility on behalf of the manufacturer shall expire.

4.4 CE Declaration of Conformity

This system was designed and manufactured based on careful selection of the harmonised standards to be observed, as well as additional technical specifications. It thus corresponds with the state-of-the-art and ensures maximum safety.

In terms of design and performance, this product complies with the European Directives and the supplementary national requirements. Conformity has been confirmed with the CE marking.

Directives:

- _ 2006/42/EG (Machinery Directive)
- _ 2011/65/EU (RoHS Directive)
- _ 2004/108/EG (EMV Directive)

Harmonised standards:

- _ DIN EN IEC 63000:2019
- _ UL 61010-1:2012
- _ CAN/CSA-C22.2 No. 61010-1:2012
- _ DIN EN 61010-1:2020
- _ DIN EN ISO 16090-1:2022
- _ ISO 8573-1

The declaration of conformity of the products is available upon request from the manufacturer.

The system has been tested according to Class A of the EMC Directive. All relevant tests were passed and confirmed by a certified inspection authority.

In general, it should be pointed out that interference emissions can occur in residential areas. This could affect:

- _ Radio services within a maximum radius of 20 m (e. g. analogue reception of radio and television, voice radio)
- _ Devices at close range (e. g. weakly shielded devices)

A health hazard can be completely excluded at any time.

4.5 Components and interfaces


Fig. 1 Overview Ceramill Matik with Ceramill Cleanstream

Ceramill Matik

- 1 Material depot
- 2 Operating element
- 3 Machining compartment

Ceramill Cleanstream

- 4 Drawer for storing accessories
- 5 Drawer for belt filter and tank
- 6 Cabinet for optional extraction system


Fig. 2 Overview of the material depot

- 1 Shelf compartments for the storage of material holders
- 2 Lighting strips for information on the status of the shelf compartments


Fig. 3 Overview of the milling compartment

- 1 Machining arm with holder clamping system

- 2 Material holder, clamped
- 3 Milling spindle
- 4 Measuring device, tool length
- 5 Tool depot, retracted
- 6 Cleaning ledge (also on opposite side)
- 7 Coolant/lubricant drain/Extraction connection of milling compartment


Fig. 4 Overview, tool depot, extended

- 1 Cover
- 2 Tool inserts
- 3 Tool holder



Fig. 5 Overview, rear side of milling system

- | | |
|--|-------------------------|
| 1 Extraction-hose connection | 14 Compressed-air inlet |
| 2 Control interface for extraction | 15 Compressed-air hose |
| 3 Key switch, emergency release of milling compartment | |
| 4 Compressed-air connection | |
| 5 Ethernet interface | |
| 6 USB interface for connecting the hand-held scanner | |
| 7 USB interface for connecting a PC | |
| 8 Blind connection for service purposes | |
| 9 Fuse compartment | |
| 10 Main switch | |
| 11 Power supply connection | |
| 12 Ventilation intake for milling compartment | |
| 13 Maintenance unit | |

-
- Diagram showing the internal components of the refrigerator, including the evaporator (1), condenser (2), and compressor (3).

Fig. 7 Overview, drawer for belt filter and tank of the Ceramill Cleanstream

- 1 Filter cartridge
- 2 Belt filter
- 3 Tank for coolant/lubricant (C/L)
- 4 Waste container

Various material holders are available for the system (see catalogue).

5 Installation

5.1 Country-specific supply voltage

NOTE:

Non-observance of the country-specific supply voltage and thus operation with the wrong fuse can cause damage to the machine.

The fuses of the Ceramill Matik and Ceramill Cleanstream are preset to a mains voltage of 230 V. When the mains voltage at the installation location is 110 - 115 V:

- ▷ On the rear of the Ceramill Matik and Ceramill Cleanstream: Pull out the fuse compartment below the mains plug (see figure 5 [9] and figure 6 [2]).
- ▷ Replace fuses T3.15 with the supplied fuses T6.3A.

5.2 Setting up

NOTE:

Possible damage of the electric lines and the C/L supply line:

- ▷ Pay attention that no cables are run over when moving or aligning the system.

- The machine is intended for use only in dry, enclosed rooms. The minimum clearances around the unit are:
 - 300 mm on both sides
 - 200 mm to the rear from the rear wall on (except for service unit). The intake/outlet of the ventilation on the rear wall may not be covered.
 - 38 mm to the top with opened material replacement door
- The minimum space requirements for the system are therefore:
 - Width: 1865 mm
 - Depth: 980 mm
 - Height: 2350 mm
- The maximum weight of the Ceramill Matik is 365 kg. Together with the Ceramill Cleanstream, the maximum weight of the system is 625 kg. The set-up surface must withstand an appropriate load.
- During operation, the room temperature must be between 18 °C and 30 °C; avoid large variations in temperature.

NOTE:

Danger of tipping over when casters are not locked!

- ▷ After installing or relocating the system, secure the machine by turning down the machine feet.

5.3 Maintenance unit and air purity during operation

The spindle is equipped with a sealing-air feature. This sealing air prevents shavings/chips and dust from entering the spindle. The service unit filters possible contamination in the sealing air and thus prevents costly machine damage.

For operation, the compressed air must fulfil the following minimum requirements:

- _ Flow rate: 400 l/min
- _ Pressure: 6 - 8 bar

Requirements on air purity:

- _ Solid contaminants:
 - Class 3; Filter grade better than 5 µm for solids
- _ Water content:
 - Class 4; Max. pressure dew point +3°C
- _ Total oil content:
 - Class 3; Max. oil content 1 mg/m³

In continuous operation, the pressure and the air-flow rate are monitored by the milling machine. When the required values are fallen below, the milling program stops. Once the correct air values are available again, the machine automatically carries out the milling program.

5.4 Requirements for the extraction system

If an extraction system not recommended by the manufacturer is used, this extraction system must fulfil the following requirements:

- _ Extraction capacity: 240 m³/h
- _ Suitability for zirconium-oxide dusts
- _ HEPA microfilter (capture 99.97 % at a particle size of 0.3 µm)

The connection of the extraction hose to the machine is made with an adapter (in delivery scope). The adapter fits extraction hoses with a 38 mm interior diameter.

WARNING:

Danger from inhalation of respirable milling waste!

- ▷ Only use tested HEPA microfilters (capture 99.97 % at a particle size of 0.3 µm).

6 Application and Operation

Instructions on application and operation are shown directly on the machine display.

WARNING:

Danger from shavings/chips or tool fragments being thrown from the machine!

- ▷ Make sure that the door of the milling machine is always closed during machining!

CAUTION:

Danger of crushing at the material replacement door!

- ▷ Use recessed grip.

CAUTION:

Danger of crushing at the drawers and door of the milling machine!

- ▷ Always open and close drawers and doors carefully.

NOTE:

Insufficient milling results during dry machining and increased wear due to milling dust!

- ▷ Do not operate the milling machine without extraction!

NOTE:

Insufficient milling/grinding results during wet machining and increased wear due to milling dust!

- ▷ Do not operate the milling machine without coolant/lubricant!

NOTE:

During dry processing, the swarf sieve prevents material residues from entering the extraction pipe and clogging the extraction system. During wet processing, shavings/chips could clog the sieve and cause coolant to overflow from the milling chamber.

- ▷ Only use the swarf sieve for dry processing.

7 Handling the Coolant/Lubricant (C/L)

This chapter is only relevant for Ceramill Matik in conjunction with Ceramill Cleanstream.

CAUTION:

Health hazards are possible through improper handling of C/L!

- ▷ Use only C/L recommended by the manufacturer!
- ▷ Observe the safety instructions and handling regulations for C/L!

C/L analysis

For legal and technical reasons, the coolant/lubricant must be analysed weekly.

- ▷ Observe country-specific laws and guidelines on coolant/lubricant!

Regular analysis of the C/L is required for proper function of the milling machine. The results of the analyses must be documented and archived.

- ▷ Check the following coolant/lubricant properties weekly:
 - Concentration
 - pH-value
 - Nitrite content

The analysis is to be carried out with test strips and refractometer. These aids are simple to handle and enable the pH-value, the nitrite content and the C/L concentration to be determined.

- ▷ Insert the reaction zone of the test strip into the C/L tank.
- ▷ Shake off excess fluid.
- ▷ After one minute, compare the reaction zone with the colour scale and read off the measured value.

The reference values and the information on measures in case of exceeding or falling below the reference values can be found on an information sheet in the delivery scope of the Ceramill Cleanstream.

8 Cleaning and Maintenance

Instructions on cleaning and maintenance are shown directly on the machine display.

NOTE:

Zirconia dust is abrasive!

- ▷ Clean the window in the door carefully, using only a soft cloth; do not scrub or scour.

NOTE:

Possible machine damage through dried-in coolant/lubricant and bound in milling dust.

After wet processing:

- ▷ Clean and dry the milling chamber prior to downtimes lasting several hours.

NOTE:

The milling chamber may only be cleaned with water.

- ▷ Do not use cleaning agents.

NOTE:

The milling compartment may only be cleaned with the specified cleaning aids (brush, cleaning cloth)!

- ▷ Do not use sponges, coarse brushes etc.

NOTE:

Danger of damaging the machine!

- ▷ Do not use compressed air, ultrasonic or a steam jet for the cleaning!

NOTE:

Only use the provided collet grease for lubrication of the collet cone. Do not use other types of grease.

9 Malfunctions, Repairs and Warranty

9.1 Malfunctions

In case of malfunctions:

- ▷ Restart the software.
- ▷ Restart the milling machine.
- ▷ Restart the PC.

9.2 Repairs

Repairs may be carried out only by trained, specialised personnel.

9.3 Warranty

The warranty complies with the statutory provisions. For more information, please refer to our General Terms of Business.

10 Environmental Protection

10.1 Packaging

In terms of packaging, the manufacturer participates in country-specific recycling systems, which ensure optimal recycling.

All packaging materials used are environmentally friendly and recyclable.

10.2 Waste equipment

Waste equipment contain materials that should be sorted for environmental-friendly recycling.

- ▷ After the service life has elapsed, dispose of the machine via the public disposal systems in an environmentally-friendly manner.

The plastic parts are labelled. Thus, they can be sorted for recycling or disposal.

10.3 Coolant/Lubricant (C/L)

- ▷ Dispose of the C/L in accordance with the applicable regulations and laws.

Proper disposal of the C/L must be documented.

11 Technical Data

Subject to changes.

Ceramill Matik

	Unit	Value
Art. No.	–	181200
Dimensions (D × W × H)	mm	780 × 1265 × 1015
Weight empty/max.:	kg	310/365
Electrical connections	V/Hz	100-230/50-60
Max. mains fluctuation	%	10
Overvoltage category	–	II
Power output	W	400
Fuse (time-delayed)	A	3.15 / 6.3
Axes	–	5
Positioning accuracy	µm	< 1
Allowable temperature range (operation)	°C	18 – 30
Max. relative humidity (non-condensing)	%	80
Ambient conditions	–	Only in enclosed rooms
Max. operating height above NHN (standard elevation zero)	m	2000
Contamination degree of the intended environment	–	2
Spindle torque	Ncm	11.2
Spindle speed	rpm	100000
Collet diameter	mm	3
Sound pressure level	dB(A)	60
Interfaces	–	USB / Ethernet / WLAN
Requirement for compressed-air connection	–	Dry, clean compressed air
_ Pressure	bar	6 – 8
_ Flow rate	l/min	400
Air purity		
_ Solid contaminants	–	Class 3; Filter grade better than 5 µm for solids
_ Water content	–	Class 4; Max. pressure dew point +3 °C
_ Total oil content	–	Class 3; max. oil content 1 mg/m ³

Tab. 1

	Unit	Value
Requirement for the extraction system	–	Suitability for zirconium-oxide dusts, HEPA microfilters (capture 99.97 % at a particle size of 0.3 µm)
_ Diameter of extraction connection	mm	38
_ Suction capacity	m ³ /h	240

Tab. 1

Ceramill Cleanstream

	Unit	Value
Art. No.	–	178660
Dimensions (D × W × H)	mm	780 × 1265 × 985
Weight empty/max.:	kg	210/260
Electrical connections	V/Hz	100-230/50-60
Max. mains fluctuation	%	10
Overvoltage category	–	II
Power output	W	200
Fuse (time-delayed)	A	3.15 / 6.3
Allowable temperature range (operation)	°C	18 - 30
Max. relative humidity (non-condensing)	%	80
Ambient conditions	–	Only in enclosed rooms
Max. operating height above NHN (standard elevation zero)	m	2000
Contamination degree of the intended environment	–	2
Interfaces	–	USB
Max. flow rate of pump	l/min	40
Max. tank capacity	l	50

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Amann Girrbach AG
Amann Girrbach GmbH