ceramill[®] motion DRS

DE	Betriebsanleitung	3 - 24
EN	User Manual	25 - 46
FR	Mode d'emploi	47 - 68
IT	Istruzioni d'uso	69 - 90
ES	Modo de empleo	91 - 111

AMANNGIRRBACH

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- Translation of the original Instruction Manual -

Table of Contents

1	Explanation of Symbols	26	8	Service (Cleaning and Maintenance) ...	40
2	General Safety Instructions	27	8.1	Weekly service	40
3	Suitable Personnel	27	8.2	Monthly service	40
4	Machine Specifications	28	8.3	Manual service	41
4.1	Delivery Scope	28	9	Settings Menu	43
4.2	Intended Use	28	9.1	Network	43
4.3	CE Declaration of Conformity	28	9.2	Locking the screen	43
4.4	Machine overview	29	9.3	Wireless network	44
4.5	Components/Interior	30	9.4	Language	44
4.6	Interfaces/Connections	31	9.5	Info	44
4.7	Operating screen	32	10	Malfunctions, Repairs and Warranty ..	45
4.8	Accessories	32	10.1	Malfunctions	45
5	Installation	33	10.2	Repairs	45
5.1	Setting up the milling machine	33	10.3	Warranty	45
5.2	Establishing the Network Connection	33	10.4	External maintenance	45
5.3	Starting-up the milling machine ...	34	11	Technical data	46
6	Application and Operation	34	12	Environmental Protection	46
6.1	Symbols on the buttons	34			
6.2	Menu overview	35			
6.3	Locking/unlocking the screen	36			
6.4	Entering text	36			
7	Operation during Production	37			
7.1	Inserting a blank	37			
7.2	Inserting/replacing a tool	37			
7.3	Changing the blank holder	38			
7.4	Transfer of milling and grinding jobs	38			
7.5	Milling zirconium oxide	39			

EN

1 Explanation of Symbols

Warning indications

Warning indications in the text are marked with a triangle and boxed.

In case of hazards through electricity, the exclamation mark in the warning triangle is substituted by a lightning bolt.

Signal words at the beginning of a warning indication specify the type and severity of the consequences, if the measures to avert the hazard are not adhered to.

- _ **NOTE** means that property damage can occur.
- _ **CAUTION** means that light to fairly serious personal injury can occur.
- _ **WARNING** means that serious personal injury can occur.
- _ **DANGER** means that life-threatening personal injury can occur.

Important information

Important information that do not lead to hazards for humans or property damage are marked with the icon on the left and are also boxed.

Other symbols on the machine

Symbol	Meaning
	USB port (socket type B)
	Network connection (Ethernet)
	Fuse

Other symbols in the Manual

Symbol	Meaning
▷	Item of an operation description
—	Item of a list
▪	Subitem of an operation description or a list
[3]	Numbers in square brackets refer to position numbers in graphics/figures

2 General Safety Instructions

When installing, starting-up and operating the machine, always observe the following safety instructions:

CAUTION:

A defective machine can lead to malfunctions!
On detection of damage or a functional defect of the machine:

- ▷ Label the machine as defective.
- ▷ Prevent further operation until the machine has been repaired.

NOTE:

Possible machine damage through insufficient cooling during milling/grinding operations!

- ▷ For milling/grinding operations, operate the machine only with sufficient quantity of coolant/lubricant and the coolant circuit switched on.

NOTE:

Damage through escaping coolant/lubricant!

- ▷ Check the tank regularly for leaks.
- ▷ Make sure that no liquid escapes from the filter cap.

NOTE:

- ▷ Switch the machine off when not in use or unsupervised for longer periods, e.g., overnight. This method of saving electrical energy also benefits the environment.

3 Suitable Personnel

NOTE:

The machine may only be unpacked, installed or relocated by trained specialist personnel, authorized by the manufacturer.

NOTE:

Starting-up and operation of the machine may only be carried out by trained specialised personnel.

EN

4 Machine Specifications

4.1 Delivery Scope

- _ Milling machine
 - _ Mains cable
 - _ Wireless barcode scanner with USB receiver
 - _ Milled-out test specimens in packaging
 - _ Network cable (5 m)
 - _ WLAN stick
 - _ USB cable
 - _ Service set for collet
 - _ Torque screwdriver for blank holder
 - _ Collet tool
 - _ 3 × 125 ml coolant/lubricant
 - _ Calibration body
 - _ Calibration pin
 - _ 4 × Filter (1 of which mounted in tank)
 - _ Tank key
 - _ Maintenance oil for axes
- After unpacking, check the machine for completeness and possible transport damages. Please claim any transport damages immediately with your supplier.

4.2 Intended Use

The milling machine is a PC-controlled milling and grinding machine for the fabrication of dentures by means of wet processing of blanks made of glass ceramics, plastics or zirconium oxide.

Using blanks and tools not approved by the manufacturer can damage the machine and make the product unusable. For such cases, the manufacturer shall assume no liability whatsoever.

Unauthorised modifications/alterations of the machine shall also void the warranty.

4.3 CE Declaration of Conformity

This product was designed and manufactured based on careful selection of the harmonised standards to be observed, as well as additional technical specifications. It thus corresponds with the state-of-the-art and ensures maximum safety. In terms of design and performance, this product complies with the European Directives and the supplementary national requirements. Conformity has been confirmed with the CE marking.

Directives:

- _ 2006/42/EG Machinery Directive
- _ 2011/65/EU RoHS II Directive
- _ 2004/108/EG EMC Directive

Harmonised standards:

- _ DIN EN IEC 63000:2019
- _ UL 61010-1:2012
- _ CAN/CSA-C22.2 No. 61010-1:2012
- _ DIN EN 61010-1:2020
- _ ISO 16090-1:2022

The product's declaration of conformity is available upon request from the manufacturer.

4.4 Machine overview


Fig. 1 Machine overview, front

- 1 Screen
- 2 Milling chamber
- 3 Tank
- 4 Storage drawer

4.5 Components/Interior


Fig. 2 Machine overview, interior

- 1 Blank holder
- 2 Workpiece
- 3 Spindle
- 4 Swarf sieve
- 5 Door
- 6 Water tank
- 7 Filter housing

4.6 Interfaces/Connections


Fig. 3 Machine overview, rear

- 1 WLAN stick (delivery scope) 
- 2 USB interfaces 
- 3 Ethernet interface 
- 4 Fuse compartment
- 5 Power supply connection
- 6 Main switch

4.7 Operating screen


Fig. 4 Start screen with icons

- 1 Restoration preview
- 2 Patient's information
- 3 Restoration information
- 4 Name of practitioner
- 5 Button *Tools*
- 6 Button *Machine*
- 7 Button *Setup*
- 8 Button *Messages*
- 9 Button *Options*
- 10 Button *Setup and Start/Stop*

The operating screen is designed as an indication display with touch-buttons. It thus enables the

- _ indication of the operating status
- _ indication of malfunctions and information
- _ operation of the milling machine

4.8 Accessories

Various holders, cutting and grinding tools are available for the milling machine.

Art. No.	Designation
For glass ceramic and composite	
179604	Cooling liquid (125 ml) + filter (10 x)
751004	Roto RFID 2.0 Diamond BLG
751006	Roto RFID 1.0 Diamond BLG
751008	Roto RFID 0.4 Diamond BLG
For zirconium oxide and plastic	
181410	Motion DRS - Zirconia milling set
751010	Roto RFID 2.5 DC BLM
751011	Roto RFID 1.0 DC BLM
751012	Roto RFID 0.6 DC BLM

Tab. 1

For milling zirconium oxide, accessory 181410 "Motion DRS - Zirconia milling set" is required.

5 Installation

5.1 Setting up the milling machine

- _ The machine is intended exclusively for use within dry, closed rooms.
- _ A clearance of at least 200 mm is to be observed sideways, to the rear and upwards. Minimum space requirements for the milling machine (incl. connections):
 - Width: 580 mm
 - Depth: 490 mm
 - Height: 580 mm
- _ The empty weight of the machine is 66 kg. The set-up surface must withstand an appropriate load.
- _ During operation, the room temperature must be between 18 °C and 30 °C; avoid large variations in temperature.
- ▷ Set up the machine on a sturdy work table or work bench (off the ground and clear of walls).
- ▷ Adjust the height of the machine by screwing the feet in or out, ensuring that the machine is positioned firmly on all four feet.
- ▷ Connect the machine to the power supply connection ([5] in Fig. 3) using the original mains cable (delivery scope).

When the machine is being connected via the USB interface ([2] in figure 3):

- ▷ Only use the USB cable from the delivery scope.
- ▷ Directly connect the cable to a USB port of the PC.

When using a USB switch/hub, the function of the machine will not be sufficiently ensured.

5.2 Establishing the Network Connection

When connecting the milling machine to the Internet, you can choose between a LAN or a WLAN connection. For the connection, the milling machine must first be initialised in the software.

When operating via a network, up to eight machines can be connected and controlled via the CAM software. The PC and the milling machine must both be in the same sub-network.

5.2.1 LAN connection

- ▷ Connect the Ethernet interface ([3] in Fig. 3) of the milling machine with the PC or with the network using a network cable.
The maximum permitted length of the network cable to the PC or to the network router is 100 m.
- ▷ Switch on the milling machine.

5.2.2 Wireless connection

A WLAN stick is included in the delivery scope of the milling machine.

- ▷ Plug the WLAN stick into a USB port on the rear of the machine ([1] in figure 3).
- ▷ Establish connection to a wireless network via the menu *Settings > Wireless network*.

5.3 Starting-up the milling machine

5.3.1 Machine

- ▷ Switch the milling machine on via the main switch.
The milling machine runs through an automatic sequence via the screen for entering the basic settings:
 - Basic settings network connection
 - Referencing axis system
 - Filling tank and flushing coolant circuit
 - Inserting swarf sieve
 - Inserting tools
- ▷ Follow the instructions on the screen.
All settings can also be entered or changed later via the menu.

6 Application and Operation

**WARNING:**

Danger from shavings/chips or tool fragments being thrown from the machine!

- ▷ Make sure that the door of the milling machine is always closed during processing!

The water for the coolant/lubricant must have drinking-water quality. The optimal hardness range is between 10 - 20 °dH (corresponds to 1.8 - 3.6 mmol total hardness per litre). If the hardness is far beyond this range, we recommend the addition of distilled water. For grinding glass-ceramic materials, it is generally recommended to use distilled water.

6.1 Symbols on the buttons

Symbol	Designation	Where?
	Tools	Start screen
	Machine	
	Settings	
	Messages	
	Opening/closing the door of the milling chamber	Varies
	Diameter, tool insert	Tool over-view
	Tool wear	

Tab. 2

6.2 Menu overview


In the start screen as the top menu level, the following submenus can be opened:

- _ Tools
 - Overview
 - Picking up a tool/Unloading a tool
- _ Machine
 - Maintenance
 - Cleaning the machine
 - Cleaning the collet
 - Cleaning the coolant tank
 - Calibrating the machine
 - Specifying the day for the weekly maintenance
 - Specifying the day for the monthly maintenance
 - Changing the blank holder
 - Moving to rest position
 - Opening the protective door
 - Service
 - Calibrating the machine (monthly)
 - Cleaning the machine (weekly)
 - Cleaning the water tank (weekly)
 - Cleaning the collet (weekly)
 - Notification history
 - List of notifications
 - Jobs in stock
 - List of unprocessed jobs
 - Job history
 - List of processed jobs
- _ Settings
 - Configuration
 - Machine
 - Barcode scanner
 - Locking the screen
 - Network
 - DHCP
 - IP address
 - Subnet mask
 - Gateway
 - Wireless network
 - List of available networks
 - Language
 - List of available languages
 - Info
 - Controller
 - HMI
 - Run times and maintenance times
- _ Notifications

6.3 Locking/unlocking the screen


To prevent unwanted actions from being carried out, the screen can be locked (e.g. during cleaning):

- ▷ Press  *Settings* button.
- ▷ Press  *Lock screen* button.


To unlock:

- ▷ Press and hold round button for approx. 3 seconds.

6.4 Entering text

If text entries are required in a menu, a keypad appears on the screen:

7 Operation during Production

7.1 Inserting a blank

- ▷ Insert the blank shown on the machine screen into the blank holder.


Fig. 5 Inserting a blank

- ▷ Insert blank into holder to the stop and hold.
- ▷ Tighten screw with the provided torque screwdriver until it “clicks”.


7.2 Inserting/replacing a tool

 At least two tools are mandatory for operation.

 The tools are colour-coded depending on application. They can be placed in any slot of the tool holder. The detection of the tool type and its wear takes place automatically via RFID.

 The tools can be inserted via context menu *Insert tool* for empty tool positions and *inserted or exchanged via Change tool* for occupied tool positions.

- ▷ In the start screen, press button  *Tools*.
The screen displays:

To insert/replace a tool:

- ▷ Select desired tool location.

The screen displays:

- ▷ Press *Take* or *Unload*.

The blank holder moves to the insertion position and the door opens.

- ▷ Insert tool into or remove from tool holder.
- ▷ Close door manually.

The automatic tool inventory is carried out.

7.3 Changing the blank holder

In delivery condition, the milling machine is equipped with one blank holder. It may be required to use a different blank holder for processing certain blocks.

To change the blank holder:

- ▷ In the start screen, press button  *Machine*.
- ▷ In the Service menu, press button  *Change blank holder* and follow the instructions on the screen.


NOTE:

Inaccurate working results or damage to the blank holder possible!

- ▷ Tighten fastening screw only with the supplied torque wrench.

7.4 Transfer of milling and grinding jobs

- ▷ Transferring jobs from the CAD software to the machine.

If the wrong job is selected, it can be deleted using the button .



- ▷ When the job is successfully transferred: Press *Setup button*.

All necessary components are checked for via the screen, e.g.:

- Is the correct blank holder inserted?
- Are the necessary tools available?
- Is the correct block inserted?
- Is the correct coolant tank inserted?

- ▷ To verify the block: Scan the material code from the block or the packaging.


- ▷ After the job is set up: Press the **START** button.
- ▷ To interrupt or stop the program: Press **STOP** button.

 At the end of the milling/grinding program, the door opens automatically.

7.5 Milling zirconium oxide

When milling zirconium oxide, an additional query is made as to whether the necessary water tank is inserted.


 To prevent contamination from one material to another, the machine distinguishes between the supplied water tank (for processing glass ceramics and plastics) and the zirconium water tank 181410.

8 Service (Cleaning and Maintenance)

The cleaning of the milling machine includes the interior, the collet and the tank.

- ▷ Carry out service as instructed by the milling machine (weekly/monthly).

A wizard on the screen of the milling machine guides you through the necessary service steps.

NOTE:

Possible damage to the milling machine!
Improper cleaning agents can damage the milling machine.

- ▷ Do not use compressed air, ultrasound or a steam jet for the cleaning (except for cleaning of the collet).
- ▷ Clean the milling chamber with water only. Do not use cleaning agents.
- ▷ Clean and dry the milling chamber prior to downtimes lasting several hours.

8.1 Weekly service

The weekly service consists of the following steps:

- _ Cleaning the interior, removing debris from the nozzles
- _ Cleaning the collet/Spindle maintenance
- _ Cleaning the tank and replacing coolant/lubricant (C/L)

8.2 Monthly service

The monthly service consists of the following steps:

- _ Cleaning the interior, removing debris from the nozzles
- _ Cleaning the collet/Spindle maintenance
- _ Cleaning the tank and replacing coolant/lubricant (C/L)
- _ Machine calibration


8.3 Manual service


Independent of the specified service intervals, manual service can be performed at any time.

► In the start screen, press button .

The screen displays:

► Press *Maintenance button* and select the desired maintenance.

As a reminder for the regular maintenance, a day of the week can be set.

8.3.1 Cleaning the machine

► Press  *Clean machine button* and follow the wizard's instructions.

8.3.2 Cleaning the collet

The provided spindle service set is required for maintenance of the spindle.


Fig. 6 Service set for collets

- 1 Collet tool
- 2 Spindle brush
- 3 Collet brush
- 4 Collet grease

NOTE:

Only use the provided collet grease for lubrication of the collet cone. Do not use other types of grease.

 For remaining debris, the collet can be cleaned in distilled water in an ultrasonic bath/with a steam jet cleaner.

► Do not use solvent for cleaning.

 To prevent contamination of the spindle, compressed air is actively expelled during collet cleaning. This may cause whistling noises.

► Press  *Clean collet button* and follow the wizard's instructions.

8.3.3 Cleaning the tank and replacing coolant/lubricant (C/L)

For cleaning the tank and replacing the coolant/lubricant, a new filter, 125 ml of coolant/lubricant and the tank key are required.

CAUTION:

Health hazards are possible through improper handling of C/L!

- ▷ Use only C/L recommended by the manufacturer!
- ▷ Observe the safety instructions and handling regulations for C/L!

Dispose of the C/L in accordance with the applicable regulations and laws.

Perform the tank cleaning over a washbasin/sink.

The water for the coolant/lubricant must have drinking-water quality. The optimal hardness range is between 10 - 20 °dH (corresponds to 1.8 - 3.6 mmol total hardness per litre). If the hardness is far beyond this range, we recommend the addition of distilled water. For grinding glass-ceramic materials, it is generally recommended to use distilled water.

- ▷ Press  *Clean coolant tank button* and follow the wizard's instructions.

8.3.4 Calibrating the machine

For calibration of the machine, the provided calibration body, the torque wrench and the calibration pin are required.

NOTE:

Faulty calibration possible!

- ▷ Clean the machine thoroughly prior to calibration (see section 8.3.1).


- ▷ Press  *Calibrate machine button* and follow the wizard's instructions.


9 Settings Menu

▷ In the start screen, press the  *Settings button*.

The screen displays the Settings menu:


9.1 Network

Network settings can be made in this menu.

▷ Press the  *Network button*.

The screen displays:


9.2 Locking the screen

To prevent unwanted actions from being carried out, the screen can be locked:


▷ Press the  *Settings button*.

▷ Press the  *Lock screen button*.

To unlock:

▷ Press and hold round button for approx. 3 seconds.

9.3 Wireless network


Connection to the Internet is described in section 5.2.

Wireless network settings can be made in this menu.

▷ Press the *Wireless network button*.

The screen displays the list of available wireless networks:

To connect to a wireless network:

- ▷ Tap on the desired network.
- ▷ Enter the network password.

9.4 Language

▷ Press the *Language button*.

The screen displays the available languages:

▷ Select the language and close the window.

The menu language is changed accordingly.

9.5 Info

▷ Press the *Info button*.

The screen displays the machine name and the screen and audio settings:

10 Malfunctions, Repairs and Warranty

10.1 Malfunctions

In case of malfunctions:

- ▷ Restart the software.
- ▷ Restart the milling machine.
- ▷ Restart the PC.

10.2 Repairs

Repairs may be carried out only by trained, specialised personnel.

10.3 Warranty

The warranty complies with the statutory provisions. For more information, please refer to our General Terms of Business.

10.4 External maintenance

Regardless of the number of operating hours, external maintenance must be carried out once a calendar year. Further information on maintenance can be obtained from your distributor or service partner.

11 Technical data

Subject to changes.

	Unit	Value
Art. No.	–	181400
Dimensions (D × W × H)	mm	490 × 580 × 580
Weight	kg	66
Electrical connections	V/Hz	100-230/50-60
Power output	W	250
Fuse 230 V / 115 V	A	T3.15 / T6.3
Axes	–	4
Accuracy	μm	≤ 5
Spindle torque	Ncm	9.2
Spindle speed	rpm	100000
Collet diameter	mm	3
Sound pressure level	dB(A)	68
Interfaces	–	USB/Ethernet
Max. pump flow rate	l/min	2.1
Ambient conditions in operation:		
– Temperature	°C	+18 ... +30
– Humidity	%	80
– Maximum operating elevation (standard elevation zero)	m	2,000

Tab. 3

	Unit	Value
Ambient conditions storage, transport:		
– Temperature	°C	-20 ... +65
– Humidity	%	80
– Ambient pressure)	mbar	500 ... 1,060

Tab. 3

12 Environmental Protection

Packaging

In terms of packaging, the manufacturer participates in country-specific recycling systems, which ensure optimal recycling.

All packaging materials used are environmentally-friendly and recyclable.

Waste equipment

Waste equipment contain materials that should be sorted for environmental-friendly recycling.

▷ After the service life has elapsed, dispose of the machine via the public disposal systems in an environmentally-friendly manner.

The plastic parts are labelled. Thus, they can be sorted for recycling or disposal.

Coolant/Lubricant (C/L)

▷ Dispose of the C/L in accordance with the applicable regulations and laws.

Proper disposal of the C/L must be documented.

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Amann Girrbach AG