

ACCESSORIES TO SUPPORT
YOUR RESCUE



WALL CABINETRY

Secure & Accessible



Keep your ZOLL AED 3 defibrillator secure and easily accessible for medical emergencies with a variety of battery-operated, alarmed wall cabinets.* Each holds a ZOLL AED 3, a backup set of electrodes, and a backup battery should you need them. Choose the surface-mounted model for easy installation or the semi-recessed or fully recessed models for less intrusive storage. All three models include installation instructions.



Fully Recessed Wall Cabinet

Part# 8000-001258

Meets ADA (Americans with Disabilities Act) requirements
14 inches x 14 inches x 1.5 inches
(35.5 cm x 35.5 cm x 3.8 cm)



Standard Surface-mounted Wall Cabinet

Part# 8000-001256

17.5 inches x 17.5 inches x 7 inches
(45.5 cm x 45.5 cm x 17.8 cm)



Semi-recessed Wall Cabinet

Part# 8000-001257

Meets ADA (Americans with Disabilities Act) requirements
14 inches x 14 inches x 3.5 inches
(35.5 cm x 35.5 cm x 8.9 cm)

*9-volt battery included

AED SIGNAGE

Create High Visibility

AED Sign Package

3D wall sign and 8.5" x 11" wall-mounted sign to improve AED visibility

Part# 8000-0825



Strobe for Standard Surface-mounted Wall Cabinet

Part# 8000-001259



Strobe for Recessed Wall Cabinet

Part# 8000-001267



WALL MOUNTS

Simple & Universal

The universal wall mount holds the ZOLL AED 3 alone or in its carry case. It includes mounting hardware, consumable tracking card, wall-mount template, and installation instructions.



Wall Mount Bracket for ZOLL AED 3 (without carry bag)

Part# 8000-001255



Wall Mount Bracket for ZOLL AED 3 with carry case and backup CPR Uni-padz electrodes (in carry case)

Part# 8000-001266

CARRY CASES

Rugged & Portable


Choose from four different portable storage options to safely carry the ZOLL AED 3 defibrillator: the standard carry case, the small hard-shell carry case, the large hard-shell carry case, or the rescue backpack. The hard-shell cases are waterproof, with an ingress protection (IP) rating of IP68 for dust and water.

Standard Carry Case

In addition to the ZOLL AED 3 with pre-connected CPR Uni-padz electrodes, the standard carry case can also hold a backup set of CPR Uni-padz. Part# 8000-001250

Replacement Shoulder Strap for ZOLL AED 3 Carry Case

Part# 8000-001252



Large Hard-shell Carry Case

Protect and carry everything needed for a vigorous rescue with the large hard-shell carry case, including a backup set of CPR Uni-padz electrodes and a backup battery. Part# 8000-001254



Small Hard-shell Carry Case

For lightweight maximum protection, choose the small hard-shell carry case. It accommodates the ZOLL AED 3 and one spare battery. Part# 8000-001253



AED Rescue Backpack

Rescue backpack for ZOLL AEDs features outer pouch that holds the ZOLL AED 3 or AED Plus®. Main compartment features three removable, stick-in-place pouches. Part# 8000-001468-01



Battery Case

Store the ZOLL AED 3 battery in this convenient pouch for easy identification and access when needed. Battery case stores easily in both the large and small hard-shell carry cases. Part# 8000-001251

ZOLL AED 3 TRAINER AND ACCESSORIES


Simulate a Rescue with Real CPR Help

The ZOLL AED 3® Trainer is the only automated external defibrillator (AED) training solution that uses Real CPR Help® to provide real-time feedback on a trainee's compression quality.

- Trainees learn the proper rate and depth of CPR chest compressions.
- An active LCD screen with full-color rescue images mirrors what a student will see on a clinical ZOLL AED 3 device.
- An on-device pause button allows the instructor to pause a training scenario.

- Trainer device ships with training harness, training electrode gels/liner, batteries, and Operator's Guide.

Part # 8028-000001-01



Replacement CPR Uni-padz Electrode Liners


Part # 8028-000013



CPR Uni-padz Electrode Training Harness

For use with ZOLL AED 3 Trainer.

Part # 8028-000010



Replacement Adhesive Gels for CPR Uni-padz Electrodes

Gel must be applied to the reusable pad before use.

Part # 8028-000014 (Case of 5)



ACCESSORIES AND SERVICES AT YOUR FINGERTIPS

Want to know more about ZOLL AED 3 accessories and services?

Contact your local authorized ZOLL distributor or visit ZOLL at [ZOLL.com/aed3](https://www.zoll.com/aed3)

Authorized ZOLL distributors can help you with your complete AED program, including:

- Program management and maintenance
- Compliance with local regulations and registries
- Training for your volunteer responders
- Replenishment of accessories
- Adding more AEDs to your program

ZOLL MEDICAL CORPORATION

269 Mill Road | Chelmsford, MA 01824 | 978-421-9655 | 800-804-4356 | [zoll.com](https://www.zoll.com)

Copyright © 2021 ZOLL Medical Corporation. All rights reserved. AED Plus, Uni-padz, Real CPR Help, ZOLL AED 3, and ZOLL are trademarks or registered trademarks of ZOLL Medical Corporation in the United States and/or other countries. All other trademarks are the property of their respective owners.

Printed in U.S.A.
MCN PP 2110 0621



[dhpsupply.com](https://www.dhpsupply.com) | 800.626.2163

ZOLL®