



AssureCheck™ 10HR


Biological Indicator

Don't wait
to validate.

In-office assurance for
same-day confidence.



dhpsupply.com
800.626.2163



Test–Incubate–Validate:
In just one day

Assured results at your fingertips—it's about time!

In-office biological monitoring for faster results and a more streamlined workflow.



Faster turnaround time

In-office biological indicator results can be available the same day, versus typical 3–7 days needed when using a mail-out option.



Immediate risk mitigation

Mail-out testing can delay response times, potentially exposing patients and staff before issues are identified, while in-house biological monitoring enables immediate responses to failed results.



Enhanced compliance

In-house testing makes it easier to adhere to testing frequency recommendations of the CDC and State Dental Boards and reduces the risk of non-compliance during audits or inspections.



Lower long-term costs

Save on shipping costs, lost time, and potential penalties for compliance lapses. While initial investment in in-house equipment may be higher, the per-test cost is typically lower over time compared to ongoing mail-out fees.

The Gold Standard in Sterilization Monitoring

Biological indicators (BIs) are the most accurate method for verifying that the conditions for sterilization have been met by confirming whether microorganisms have been killed during a sterilization cycle. For dental practices, this is critical in ensuring the safety of your patients and practice.



AssureCheck™ 10HR
ST1010, 25 units/box; 6 boxes/case

AssureCheck™ 10HR biological indicators deliver easy-to-read results the same day, so you can maintain compliance and protect your patients with confidence.

- **10-hour final result**
Be assured within a standard workday.
- **Clear visual readout**
Easy-to-see color change from purple to yellow.
- **Versatile**
Compatible for use in all steam sterilizers.



AssureCheck™ Dry Block Incubator is lightweight, compact, and easy to use. Pre-programmed for optimal temperature and designed to work seamlessly with AssureCheck™ 10HR biological indicators, this incubator lets you take control of your sterilization monitoring process.



AssureCheck™ Dry Block Incubator
ST1000, 1/box

W 3.5" / 9 cm D 4.0" / 10.2 cm
H 2.0" / 5.08 cm W 0.89 lb / 0.4 kg

Pre-set temperature
No calibration needed

Built in vial crusher


Dual well design
Accommodates control and test BIs

Take control of your sterilization monitoring process and eliminate the dependency on external laboratories.
With AssureCheck™ 10HR biological indicators and dry block incubator – you can.



- Streamline compliance
- Fast-track your results—
10 hours instead of days
- Empower your dental team

Ready to take on in-office sterilization testing?
Scan for more info and helpful tools.



dhpsupply.com | 800.626.2163

medicom.com | 1-800-361-2862

Rev.00
©2025 AMD Medicom Inc. All rights reserved.

