

Design meets Efficiency



VITA SMART.FIRE® ADVANCED

The compact and modern
furnace for dental practices
and laboratories with space
for bridges with up to 14 units.



For all common materials.

VITA SMART.FIRE® ADVANCED furnace

With firing temperatures of up to 1200 °C, the VITA SMART.FIRE ADVANCED is suitable for

- veneering, crystallization, glaze-, stain- and corrective firings
- crystallization firings with and without vacuum

Thanks to the preprogrammed firing programs, which are optimised for VITA ceramic materials, it not only offers efficient processing, but also esthetically pleasing results.

