

STERILIZATION AND DRYNESS

at speeds faster than conventional chambered autoclaves



001_ 07

East Cycle Times

STATIM 2000 G4 unwrapped cycle times are the world's fastest at 6 minutes. Wrapped cycles, including drying, can be completed in 35 minutes.

STATIM 5000 G4 can complete unwrapped cycles in as little as 9 minutes. Wrapped cycles, including drying, can be completed in 40 minutes.

G4 Technology

SciCan's G4 equipment can connect and send information to your smart devices so you can have cycle data, instructions, and maintenance notifications at your fingertips.



Reliable

For over 25 years, STATIM's unmatched speed, reliability, and dependability have made it one of the world's best selling cassette autoclaves.



Capacity

STATIM G4 autoclaves not only offer fast sterilization times but also enough capacity to meet the needs of a busy practice.

STATIM 2000 G4 capacity allows for up to 5 pouches and the STATIM 5000 G4 allows for up to 12 pouches.



Gentle

The fresh water feature guarantees fresh water with every cycle. State of the art steam technology is gentle on both solid and hollow instruments.

STATIM 2000 G4 Specifications

p/n G4-121101							
Unit size (L x W x H): 19.5" x 16.25" x 6" / 49.5 cm x 41.5 cm x 15 cm Cassette internal dimensions: 11" x 7" x 1.4" / 28 cm x 18 cm x 3.5 cm Minimum shelf height for proper leveling: 8.5" / 22 cm Weight without water: 46 lbs / 21 kg		Unwrapped	Wrapped	Rubber / Plastic	Air Drying		
	Sterilization Time	3.5 minutes	10 minutes	15 minutes	Air drying only		
	Sterilization Temp	135° C	135° C	121º C	Air drying only		
	Total Cycle Time*	6 minutes	14 minutes	20 minutes	Air drying only		
Power consumption:	*Ousis time is dependent an instrument lead and does not include during						

*Cycle time is dependent on instrument load and does not include drying.

STATIM 5000 G4 Specifications

p/11 G4=201103							
Unit size (L x W x H): $23.6" \times 16.1" \times 7.5" / 60 \text{ cm} \times 41 \text{ cm} \times 19 \text{ cm}$ Cassette internal dimensions: $15" \times 7" \times 3" / 38 \text{ cm} \times 18 \text{ cm} \times 7.5 \text{ cm}$ Minimum shelf height for proper leveling: 10" / 26 cm Weight without water: 73 lbs / 33 kg		Unwrapped	Wrapped	Rubber / Plastic	Heavy Duty Unwrapped		
	Sterilization Time	3.5 minutes	6 minutes	35 minutes	6 minutes		
	Sterilization Temp	132° C	132° C	121° C	132° C		
	Total Cycle Time*	9 minutes	17.5 minutes	45 minutes	12 minutes		
Power consumption:	*Cycle time is dependent on instrument load and does not include drying.						

115 V, 60 Hz, 11A

115 V, 60 Hz, 11A



VistaPure[™] Water Filtration

Provides two grades of high-quality water for autoclaves, dental bottles, ultrasonic cleaners and the rinse cycle of instrument washers. V3000



VistaPure[™] AutoFill

Replenish any new or older-model STATIM and BRAVO autoclave with distilled-quality water automatically, enabling truly hands-free convenience.

S3020 - AutoFill for STATIM S3030 - AutoFill for BRAVO

STATMAT[™]

silicone mat.

Protect your STATIM

01-112867S - For STATIM 2000

01-112868S - For STATIM 5000

Securely seat cassettes atop your

STATIM with a safe and hygienic



ENSURE[™] Monitoring for steam sterilization

An easy to use, in-office monitoring system for steam sterilization.

SCI-BISK-120 - Ensure BI Starter Kit (25 BI's, Incubator & Log Book)



VistaCool™ Direct-to-Drain System

Protects cabinets and drains-and saves valuable staff time-by cooling autoclave wastewater and sending it directly down the drain automatically.

V7501 - For One Autoclave V7502 - For Two Autoclaves

Manufactured by:

SciCan Ltd. 1440 Don Mills Rd., Toronto, ON, Canada, M3B 3P9 Tel: 416-445-1600 Toll free: 1-800-667-7733 customerservice@scican.com

Distributed by:

SciCan Inc. 701 Technology Dr., Canonsburg, PA 15317, USA Tel: 724-820-1600 Toll Free: 1-800-572-1211 uscustomerservice@scican.com

VistaPure[™], VistaCool[™] and DrainHub[™] are manufactured by:

Vista Research Group PO Box 321, Ashland, Ohio 44805-0321 USA Tel: (419) 285-3927 VistaResearchGroup.com

f 🗾 🛗

